



Germ TerminatorTM PRO

SALON IMPLEMENT SANITIZER

THE ONLY PRODUCT OF ITS KIND.

Revolutionary device kills bacteria, viruses, fungi, molds and yeasts.

Worldwide Patents & Patents Pending

- Simple to Use
- Takes Only Minutes
- Kill Germs (99.999%)
Using Only Water!
- Uses "STEAM" and
"DRY HEAT" to Sanitize
- Clean Chamber Provides
Hygienic Storage
Environment

SALON ASEPTICTM
CONTROL SYSTEM



BACKSCRATCHERS[®]
SALON SYSTEMS, INC.



Germ TerminatorTM

STEAM DRY

Takes Only 3 Easy Steps.



1 Pre-clean tools, remove cover and place tools in basket.



2 Pour water into cup (distilled water is recommended).



3 Replace cover and pour water into fill chamber.

A Revolutionary Breakthrough... The Germ Terminator™ PRO. Ultimate Clean

The Germ Terminator™ PRO is an ingenious and incredibly easy-to-use device designed to protect your long-term health by killing germs that breed on your salon implements. The Germ Terminator™ uses an exclusive, patented process that destroys germs through the scientifically effective and time-tested use of steam, combined with a proprietary automatic drying process. This groundbreaking device kills germs using steam (212°F, 100°C) and stores implements in the Germ Terminator™ clean, dry heat compartment. The sanitized tools stand ready for your next customer. All you need to do is fill the Germ Terminator™ water reservoir with distilled water. Filling the reservoir automatically turns on the Germ Terminator™ “steam” mode. The “dry” mode is activated automatically as soon as the “steam” cycle is complete.



The Germ Terminator™ PRO Clean Chamber—The Safest Place to Store Your Salon Tools

Another important protective feature of the Germ Terminator™ is its dry-heated clean chamber, that guards each implement against cross-microbial contamination as long as the chamber has not been opened after the “dry” mode is complete. Once an implement has been sanitized in the Germ Terminator™, it remains sanitized in the clean chamber and ready for next use.

Today’s savvy consumers are becoming more knowledgeable and concerned about their own health within the salon. Use your Germ Terminator™ for every customer, every day. How many products do you own that get that kind of use? Sanitize your implements—it just makes sense.



Takes Only Seconds to Use

The Germ Terminator™ PRO is the most effective, easy-to-use device ever created to sanitize your salon implements after every use! It takes only seconds of your time per use. First pre-clean your implements, then place them in the Germ Terminator™, position the top cover, and add distilled water—it's that simple. The Germ Terminator™ uses time-tested steam and dry heat for consistently thorough sanitization.

With the Germ Terminator™ there is no exposure to other harmful disinfecting chemicals. The Germ Terminator™ is not only the most effective method of sanitization (with a kill rate of greater than 99.9999% on all germs tested), it is also the safest. Without proper implement hygiene, you may re-infect yourself, and your customers.

Independent Laboratory Testing

Independent labs performed microbiological testing on the effect that the Germ Terminator™ had with the following representative species of bacteria, viruses and fungi, including molds and yeasts:

1. Human Immunodeficiency Virus (HIV) - Type 1
2. Herpes Simplex Virus - Type 1
3. Hepatitis C Virus
4. Candida albicans
5. Staphylococcus aureus
6. Coli
7. Salmonella choleraesuis
8. Listeria monocytogenes
9. Porphyromonas gingivalis
10. Neisseria meningitidis
11. Influenza A Virus (Common Flu Virus)
12. Rhinovirus (Common Cold Virus)

List is of the most common species.

Test results from the laboratories were submitted to the United States Food and Drug Administration as part of the Germ Terminator™ FDA 510(k) medical device application.

SALON ASEPTIC™

CONTROL SYSTEM

ANOTHER FINE PRODUCT FROM:
BACKSCRATCHERS SALON SYSTEMS, INC.
9068 Elkmont Way, Elk Grove, CA 95624 USA
www.backscratchers.com

OTHER SALON SANITATION PRODUCTS AVAILABLE:
Surface cleaner, liquid disinfectant, waterless skin cleanser, soaking and storage trays, salon placemats and towels

The GERM TERMINATOR™ PRO
www.backscratchers.com



*The Germ Terminator™ PRO is a United States of America FDA-regulated medical device.

Medical device manufacturers are required to submit a 510(k) to the United States Food & Drug Administration and receive subsequent clearance from the FDA if they intend to introduce a medical device into commercial distribution. (See the FDA website at www.fda.gov for more information.)

FDA 510(k)-Cleared Medical Device
(For use worldwide)
Intl: 110v 60hz-220